

Anti-Transferrin Receptor/TFRC Antibody Picoband® Biotin Conjugated

Catalog Number: PB9233-Biotin

About TFRC

Transferrin receptor protein 1 (TfR1), also known as Cluster of Differentiation 71 (CD71), is a protein that in humans is encoded by the TFRC gene. It is mapped to 3q29. TFRC is a transmembrane glycoprotein composed of two disulfide-linked monomers joined by two disulfide bonds. Expression of human TFR1 in hamster cell lines markedly enhanced the infection of viruses pseudotyped with the glycoprotein of Machupo, Guanarito, and Junin viruses. TFR1 is a cellular receptor for New World hemorrhagic fever arenaviruses. It is required for iron delivery from transferrin to cells.

Overview

Product Name	Anti-Transferrin Receptor/TFRC Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P02786

Technical Details

Immunogen	E.coli-derived human TFRC recombinant protein (Position: M1-N198). Human TFRC shares 72% and 70% amino acid (aa) sequences identity with mouse and rat TFRC, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

2 Publications Citing This Product

1. PubMed ID: 10.3892/or.2019.7189, Sulfasalazine-induced ferroptosis in breast cancer cells is reduced by the inhibitory effect of estrogen receptor on the transferrin receptor

2. PubMed ID: 28008944, Amyloid precursor protein modulates Nav16 sodium channel currents through a Go-coupled JNK pathway

Visit bosterbio.com/anti-tfrc-picoband-trade-antibody-pb9233-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transferrin Receptor/TFRC Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.