

## Anti-p73/TP73 Antibody Picoband® FITC Conjugated

Catalog Number: PB9235-FITC

### About TP73

p73, also known as tumor protein 73 (TP73), encodes a member of the p53 family of transcription factors involved in cellular responses to stress and development. It is mapped to a region on chromosome 1p36 that is frequently deleted in neuroblastoma and other tumors, and thought to contain multiple tumor suppressor genes. The demonstration that this gene is monoallelically expressed (likely from the maternal allele), supports the notion that it is a candidate gene for neuroblastoma. Furthermore, recent findings are suggesting that over-expression of transcription factors involved in cell cycle regulation and synthesis of DNA in mammalian cells (e.g.: E2F-1) induces the expression of p73. In addition, p73 is a substrate of the c-Abl kinase and that the ability of c-Abl to phosphorylate p73 is markedly increased by gamma-irradiation.

### Overview

Product Name	Anti-p73/TP73 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O15350

### Technical Details

Immunogen	E.coli-derived human p73 recombinant protein (Position: M1-E198). Human p73 shares 96% amino acid (aa) sequence identity with mouse p73.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: 15884108, Expression and clinical significance of TAp73?, p53, PCNA and apoptosis in hepatocellular carcinoma
2. PubMed ID: 20880063, Guo H, Nan K, Hu T, Meng J, Hui W, Zhang X, Qin H, Sui C. Hepatol Res. 2010 Nov;40(11):1107-16. Doi: 10.1111/J.1872-034X.2010.00721.X. Epub 2010 Sep 28. Prognostic Significance Of Co-Expression Of Nm23 And P57 Protein In Hepatocellular Carcinoma.

Visit [bosterbio.com/anti-p73-picoband-trade-antibody-pb9235-boster.html](http://bosterbio.com/anti-p73-picoband-trade-antibody-pb9235-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p73/TP73 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.