

## Anti-TNF alpha Antibody Picoband® Cy3 Conjugated

Catalog Number: PB9246-Cy3

### About Tnf

TNFalpha (Tumor Necrosis Factor alpha) gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR2. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, and cancer. Knockout studies in mice also suggested the neuroprotective function of this cytokine.

### Overview

Product Name	Anti-TNF alpha Antibody Picoband® Cy3 Conjugated
Reactive Species	Mouse
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P06804

### Technical Details

Immunogen	E.coli-derived mouse TNF alpha recombinant protein (Position: L80-Q227). Mouse TNF alpha shares 79% and 95% amino acid (aa) sequence identity with human and rat TNF alpha, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 71 Publications Citing This Product

1. PubMed ID: 33607154, Liao L,Huang L,Wei X,Yin L,Wei X,Li T.Bioinformatic and biochemical studies of formononetin against liver injure.Life Sci.2021 Feb 16;119229.doi:10.1016/j.lfs.2021.119229.Epub ahead of print.PMID:33607154.

2. PubMed ID: 33607154, Liao L,Huang L,Wei X,Yin L,Wei X,Li T.Bioinformatic and biochemical studies of formononetin against liver injure.Life Sci.2021 Feb 16;119229.doi:10.1016/j.lfs.2021.119229.Epub ahead of print.PMID:33607154.

3. PubMed ID: 33642941, Ye D,Hu Y,Zhu N,Gu W,Long G,Tao E,Fang M,Jiang M. Exploratory Investigation of Intestinal Structure and Function after Stroke in Mice. Mediators Inflamm.2021 Feb 15;2021:1315797.doi: 10.1155/2021/1315797.PMID:33642941;PMCID:PMC7902147.

Visit [bosterbio.com/anti-tnf-alpha-picoband-trade-antibody-pb9246-boster.html](https://bosterbio.com/anti-tnf-alpha-picoband-trade-antibody-pb9246-boster.html) to see all 71 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNF alpha Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.