

## Anti-Fas Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9252-Fluoro488

### About FAS

FAS (also known as surface antigen APO1 or CD95) is a member of the tumour-necrosis receptor factor family of death receptors. It acts as an inducer of both neurite growth in vitro and accelerated recovery after nerve injury in vivo. FAS antigen is expressed and functional on papillary thyroid cancer cells and this may have potential therapeutic significance. The FAS antigen shows structural homology with a number of cell surface receptors, including tumor necrosis factor (TNF) receptors and the low-affinity nerve growth factor receptor (NGFR) and it is mapped to 10q24.1. The FAS and FASL system plays a key role in regulating apoptotic cell death and corruption of this signalling pathway has been shown to participate in immune escape and tumorigenesis.

### Overview

Product Name	Anti-Fas Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P25445

### Technical Details

Immunogen	E.coli-derived human Fas recombinant protein (Position: Q26-N173). Human Fas shares 55% and 59% amino acid (aa) sequences identity with mouse and rat Fas, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 24 Publications Citing This Product

1. PubMed ID: PMID:26617948, Significance of Fas and FasL protein expression in cardiac carcinoma and local lymph node tissues
2. PubMed ID: 10.3748/wjg.v16.i10.1274, Clinical significance of Fas and FasL protein expression in gastric carcinoma and local lymph node tissues
3. PubMed ID: 10.1039/C6FO00071A, The effect of geniste on Abeta25-35-induced PC12 cell apoptosis through the JNK-dependent Fas pathway

Visit [bosterbio.com/anti-fas-picoband-trade-antibody-pb9252-boster.html](http://bosterbio.com/anti-fas-picoband-trade-antibody-pb9252-boster.html) to see all 24 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fas Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.