

## Anti-Integrin alpha 5/ITGA5 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9254-Fluoro647

### About ITGA5

Integrin alpha-5, also known as FNRA or VLA5A, is a protein that in humans is encoded by the ITGA5 gene. The product of this gene belongs to the integrin alpha chain family. Integrins are integral membrane proteins composed of an alpha chain and a beta chain. This gene encodes the integrin alpha 5 chain. Alpha chain 5 undergoes post-translational cleavage in the extracellular domain to yield disulfide-linked light and heavy chains that join with beta 1 to form a fibronectin receptor. In addition to adhesion, integrins are known to participate in cell-surface mediated signalling.

### Overview

Product Name	Anti-Integrin alpha 5/ITGA5 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P08648

### Technical Details

Immunogen	E.coli-derived human Integrin alpha 5 recombinant protein (Position: E787-E992). Human Integrin alpha 5 shares 86% amino acid (aa) sequence identity with mouse Integrin alpha 5.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

---

1. PubMed ID: 19918321, Characterization of Rat Hair Follicle Stem Cells Selected by Vario Magnetic Activated Cell Sorting System
2. PubMed ID: 26680220, Silencing of Paternally Expressed Gene 10 Inhibits Trophoblast Proliferation and Invasion
3. PubMed ID: 24643246, A Small Physiological Electric Field Mediated Responses of Extravillous Trophoblasts Derived from HTR8/SVneo Cells: Involvement of Activation of Focal Adhesion Kinase Signaling

Visit [bosterbio.com/anti-integrin-alpha-5-picoband-trade-antibody-pb9254-boster.html](http://bosterbio.com/anti-integrin-alpha-5-picoband-trade-antibody-pb9254-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alpha 5/ITGA5 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.