

## Anti-Kallikrein 4/KLK4 Antibody Picoband® Cy3 Conjugated

Catalog Number: PB9260-Cy3

### About KLK4

Kallikrein-related peptidase 4, also known as PSTS or EMSP1, is a protein which in humans is encoded by the KLK4 gene. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. In particular, they may serve as biomarkers for both prostate cancer and breast cancer. It has been found that KLK4 is a proliferative factor with effects on gene expression and it may have a role in prostate cancer development and progression.

### Overview

Product Name	Anti-Kallikrein 4/KLK4 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y5K2

### Technical Details

Immunogen	E.coli-derived human Kallikrein 4 recombinant protein (Position: I31-S254).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Kallikrein 4/KLK4 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.