

Anti-Mitofusin 2/MFN2 Antibody Picoband® APC Conjugated

Catalog Number: PB9265-APC

About MFN2

Mitofusin-2 is a protein that in humans is encoded by the MFN2 gene. It is mapped to chromosome 1 and encodes a 757-amino acid protein that contains an ATP/GTP-binding site motif. This gene is expressed in many tissues and cell lines such as brain and KG-1 with the highest expression in heart and skeletal muscle. It has been found that MFN2 triggers mitochondrial energization, at least in part, by regulating OXPHOS expression through signals that are independent of its role as a mitochondrial fusion protein. And it contributes to the maintenance and operation of the mitochondrial network. Axonal CMT type 2A and autosomal dominant HMSN VI are caused by MFN2 and mutations in MFN2, which emphasizes its important role of mitochondrial function for both optic atrophies and peripheral neuropathies.

Overview

Product Name	Anti-Mitofusin 2/MFN2 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O95140

Technical Details

Immunogen	E.coli-derived human Mitofusin 2 recombinant protein (Position: V601-R757). Human Mitofusin 2 shares 96% and 95% amino acid (aa) sequence identity with mouse and rat Mitofusin 2, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

2 Publications Citing This Product

1. PubMed ID: 32884840, Jiao Z,Wu Y,Qu S.Fenpropathrin induces degeneration of dopaminergic neurons via disruption of the mitochondrial quality control system.Cell Death Discov.2020 Aug 25;6:78.doi:10.1038/s41420-020-00313-y.PMID:32884840;PMCID:PMC7447795.

2. PubMed ID: 24928681, Montaigne D, Marechal X, Coisne A, Debry N, Modine T, Fayad G, Potelle C, El Arid Jm, Mouton S, Sebti Y, Duez H, Preau S, Remy-Jouet I, Zerimech F, Koussa M, Richard V, Nevieri R, Edme JI, Lefebvre P, Staels B. Circulation. 2014 Aug 12;130(7):554-...

Visit bosterbio.com/anti-mitofusin-2-picoband-trade-antibody-pb9265-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mitofusin 2/MFN2 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.