

Anti-LKB1/STK11 Antibody Picoband® APC Conjugated

Catalog Number: PB9270-APC

About STK11

Serine/threonine kinase 11 or LKB1 is a protein kinase which in humans is encoded by the STK11 gene. The STK11 gene, which encodes a member of the serine/threonine kinase, regulates cell polarity and functions as a tumour suppressor. STK11 is a primary upstream kinase of adenine monophosphate-activated protein kinase (AMPK), a necessary element in cell metabolism that is required for maintaining energy homeostasis. It has been found that LKB1 physically associates with p53 and regulates specific p53-dependent apoptosis pathways. The STK11 gene extends over 23 kb of genomic DNA and is composed of 9 exons, which are transcribed in telomere-to-centromere direction.

Overview

Product Name	Anti-LKB1/STK11 Antibody Picoband® APC Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q15831

Technical Details

Immunogen	E.coli-derived human LKB1 recombinant protein (Position: K62-C430). Human LKB1 shares 89% amino acid (aa) sequence identity with both mouse and rat LKB1.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 25973293, Tang Y, Li M, Wang YI, Threadgill Md, Xiao M, Mou Cf, Song Gl, Kuang J, Yang X, Yang L, Gao Xj, Wang Yp, Meng Yp. Am J Cancer Res. 2015 Jan 15;5(2):498-513. Ecollection 2015. Art1 Promotes Starvation-Induced Autophagy: A Possible Protective Role I...

Visit bosterbio.com/anti-lkb1-picoband-trade-antibody-pb9270-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LKB1/STK11 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.