

Anti-Von Willebrand Factor/VWF Antibody Picoband® PE Conjugated

Catalog Number: PB9273-PE

About VWF

Von Willebrand factor (VWF) is a blood glycoprotein involved in hemostasis. It is mapped to 12p13.31. The VWF gene encodes von Willebrand factor (VWF), a large multimeric glycoprotein that plays a central role in the blood coagulation system, serving both as a major mediator of platelet-vessel wall interaction and platelet adhesion, and as a carrier for coagulation factor VIII. VWF released from endothelial cell Weibel-Palade bodies bound particularly avidly to the extracellular matrix. VWF deficiency or dysfunction (von Willebrand disease) leads to a bleeding tendency, which is most apparent in tissues having high blood flow shear in narrow vessels.

Overview

Product Name	Anti-Von Willebrand Factor/VWF Antibody Picoband® PE Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8CIZ8

Technical Details

Immunogen	E.coli-derived mouse VWF recombinant protein (Position: M1304-E1452).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

37 Publications Citing This Product

1. PubMed ID: PMID:26191130, Calcium signaling is implicated in the diffuse axonal injury of brain stem
2. PubMed ID: 10.12659/MSM.915442, Identification and Analysis of Differentially Expressed Genes in Human Saphenous Vein Endothelial Cells Overexpressing Domain-Containing mTOR-Interacting Protein (DEPTOR) by RNA-Seq
3. PubMed ID: 10.3892/etm.2014.1774, Construction of tissue-engineered bone using a bioreactor and platelet-rich plasma

Visit bosterbio.com/anti-vwf-picoband-trade-antibody-pb9273-boster.html to see all 37 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Von Willebrand Factor/VWF Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.