

## Anti-Lamin A+C/LMNA Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9280-Fluoro647

### About LMNA

Lamins are structural protein components of the nuclear lamina, a protein network underlying the inner nuclear membrane that determines nuclear shape and size. There are three types of lamins, A, B and C. The lamin A/C (LMNA) gene contains 12 exons. Alternative splicing within exon 10 gives rise to two different mRNAs that code for pre-lamin A and lamin C. Lamin A/C is mapped to 1q21.2-q21.3 and mutations in this gene cause a variety of human diseases including Emery-Dreifuss muscular dystrophy, dilated cardiomyopathy, and Hutchinson-Gilford progeria syndrome. Lamin A/C deficiency is thus associated with both defective nuclear mechanics and impaired mechanically activated gene transcription.

### Overview

Product Name	Anti-Lamin A+C/LMNA Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P02545

### Technical Details

Immunogen	E.coli-derived human Lamin A/C recombinant protein (Position: Y481-Y646). Human Lamin A/C shares 90% and 92% amino acid (aa) sequence identity with mouse and rat Lamin A/C, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 5 Publications Citing This Product

---

1. PubMed ID: 10.1016/j.prp.2014.11.008, Loss of lamin A but not lamin C expression in epithelial ovarian cancer cells is associated with metastasis and poor prognosis
2. PubMed ID: 33859765, Li D,Liu J,Yang C,Tian Y,Yin C,Hu L,Chen Z,Zhao F,Zhang R,Lu A,Zhang G,Qian A.Targeting long noncoding RNA PMIF facilitates osteoprogenitor cells migrating to bone formation surface to promote bone formation during aging.Theranostics.2021 Mar 20;11(11):5585-5604.doi:10.7150/thno.54477. PMID:33859765;PMCID:PMC8039942.
3. PubMed ID: 25240403, Liu J, Ma J, Wu Z, Li W, Zhang D, Han L, Wang F, Reindl Km, Wu E, Ma Q. BMC Cancer. 2014 Sep 20;14:686. Doi: 10.1186/1471-2407-14-686. Arginine Deiminase Augments The Chemosensitivity Of Argininosuccinate Synthetase-Deficient Pancreatic Cancer Cel...

Visit [bosterbio.com/anti-lamin-a-c-picoband-trade-antibody-pb9280-boster.html](http://bosterbio.com/anti-lamin-a-c-picoband-trade-antibody-pb9280-boster.html) to see all 5 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Lamin A+C/LMNA Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.