

Anti-MMP14 Antibody Picoband®

Catalog Number: PB9283

About MMP14

Matrix metalloproteinase-14, also called MT-MMP 1, is an enzyme that in humans is encoded by the MMP14 gene. The protein encoded by this gene is a member of the membrane-type MMP subfamily. This gene is mapped to 14q11.2. Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. This gene may be involved in actin cytoskeleton reorganization by cleaving PTK7. It acts as a positive regulator of cell growth and migration via activation of MMP15.

Overview


Product Name	Anti-MMP14 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-MMP14 Antibody Picoband® catalog # PB9283. Tested in WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , and 0.05 mg Na ₃ N. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50281

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human MMP14, different from the related mouse and rat sequences by one amino acid.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.

Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human

Anti-MMP14 Antibody Picoband® (PB9283) Images



Anti-MMP14 Picoband antibody, PB9283, Western blottingAll lanes: Anti MMP14 (PB9283) at 0.5ug/mlLane 1: Human Placenta Tissue Lysate at 50ugLane 2: HELA Whole Cell Lysate at 40ugLane 3: U87 Whole Cell Lysate at 40ugPredicted bind size: 55KDObserved bind size: 55KD

9 Publications Citing This Product

1. PubMed ID: 27453531, The inhibition of calpains ameliorates vascular restenosis through MMP2/TGF- β 1 pathway
2. PubMed ID: 24573354, Xiao H, Li H, Yu G, Xiao W, Hu J, Tang K, Zeng J, He W, Zeng G, Ye Z, Xu H. Oncol Rep. 2014 Apr;31(4):1832-8. Doi: 10.3892/Or.2014.3048. Epub 2014 Feb 24. Microrna-10B Promotes Migration And Invasion Through Klf4 And Hoxd10 In Human Bladder Cancer.
3. PubMed ID: 25673203, Shi Ch, Wang Ws, Zhang Zd, Li Cj, Guo Fj, Li F, Chen Am. J Huazhong Univ Sci Technolog Med Sci. 2015 Feb;35(1):111-6. Doi: 10.1007/S11596-015-1398-1. Epub 2015 Feb 12. Effect Of Lentivirus-Mediated Upa Silencing On The Proliferation And Apoptosis ...

Visit bosterbio.com/anti-mmp14-picoband-trade-antibody-pb9283-boster.html to see all 9 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MMP14 Antibody

For Research Use Only. Not for use in diagnostic procedures.