

## Anti-Nkx2.5/NKX2-5 Antibody Picoband® FITC Conjugated

Catalog Number: PB9293-FITC

### About NKX2-5

Homeobox protein Nkx-2.5, also known as NKX2E or CSX is a protein that in humans is encoded by the NKX2-5 gene. It is mapped to 5q35.1. Homeobox-containing genes play critical roles in regulating tissue-specific gene expression essential for tissue differentiation, as well as determining the temporal and spatial patterns of development. Nkx2.5 and Tbx5 directly bound to the promoter of the gene encoding cardiac-specific natriuretic peptide precursor type A (NPPA) in tandem, and both transcription factors showed synergistic activation. The cardiac homeobox protein Nkx2.5 is essential in cardiac development, and mutations in CSX (which encodes Nkx2.5) cause various congenital heart malformations.

### Overview

Product Name	Anti-Nkx2.5/NKX2-5 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P52952

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Nkx2.5, different from the related mouse and rat sequences by five amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nkx2.5/NKX2-5 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.