

## Anti-Neuropeptide Y/NPY Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9296-Fluoro550

### About NPY

This gene encodes a neuropeptide that is widely expressed in the central nervous system and influences many physiological processes, including cortical excitability, stress response, food intake, circadian rhythms, and cardiovascular function. The neuropeptide functions through G protein-coupled receptors to inhibit adenylyl cyclase, activate mitogen-activated protein kinase (MAPK), regulate intracellular calcium levels, and activate potassium channels. A polymorphism in this gene resulting in a change of leucine 7 to proline in the signal peptide is associated with elevated cholesterol levels, higher alcohol consumption, and may be a risk factor for various metabolic and cardiovascular diseases. The protein also exhibits antimicrobial activity against bacteria and fungi.

### Overview

Product Name	Anti-Neuropeptide Y/NPY Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01303

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Neuropeptide Y, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 8 Publications Citing This Product

1. PubMed ID: 10.1016/j.ijcard.2004.12.065, Expression changes of thrombospondin-1 and neuropeptide Y in myocardium of STZ-induced rats
2. PubMed ID: 10.1016/j.ijcard.2017.08.011, Increased inflammation promotes ventricular arrhythmia through aggravating left stellate ganglion remodeling in a canine ischemia model
3. PubMed ID: 10.1111/are.12979, An immunohistochemical study on endocrine cells in the neuroendocrine system of the digestive tract of milkfish *Chanos chanos* (Forsskal, 1775)

Visit [bosterbio.com/anti-neuropeptide-y-picoband-trade-antibody-pb9296-boster.html](http://bosterbio.com/anti-neuropeptide-y-picoband-trade-antibody-pb9296-boster.html) to see all 8 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Neuropeptide Y/NPY Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.