

## Anti-P2X5/P2RX5 Antibody Picoband® Biotin Conjugated

Catalog Number: PB9305-Biotin

### About P2RX5

P2X purinoceptor 5 is a protein that in humans is encoded by the P2RX5 gene. The product of this gene belongs to the family of purinoceptors for ATP. It has been identified that the P2RX5 gene within a BAC contig covering chromosome 17p13. This receptor functions as a ligand-gated ion channel. Several characteristic motifs of ATP-gated channels are present in its primary structure, but, unlike other members of the purinoceptors family, this receptor has only a single transmembrane domain. Quantitative dot blots revealed that P2X5R expression is highest in brain and immune system and is developmentally regulated.

### Overview

Product Name	Anti-P2X5/P2RX5 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q93086

### Technical Details

Immunogen	E.coli-derived human P2X5 recombinant protein (Position: D333-T422). Human P2X5 shares 42% amino acid (aa) sequence identity with rat P2X5.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.

ELISA, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 25561285, Meng M, Zheng J, Yan J, Li Q, Fang Q, Li W. Cell Biochem Biophys. 2015 Jan 6. [Epub Ahead Of Print] P2X2 And P2X5 Receptors Mediate Bladder Hyperesthesia In Icc In Female Overactive Bladder.

Visit [bosterbio.com/anti-p2x5-picoband-trade-antibody-pb9305-boster.html](http://bosterbio.com/anti-p2x5-picoband-trade-antibody-pb9305-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-P2X5/P2RX5 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.