

Anti-Parkin/PRKN Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9307-Fluoro647

About PRKN

Parkin is a RING domain-containing E3 ubiquitin ligase involved in proteasome-dependent degradation of proteins. It is mapped to 6q26. This gene is important for mitochondrial quality control by lysosome-dependent degradation of damaged mitochondria through autophagy, or mitophagy. Parkin is expressed in neuronal processes and cell bodies of neurons, but not glial cells, in the midbrain, basal ganglia, cerebral cortex, and cerebellum. Parkin associated with actin filaments, suggesting that it is a cytoskeletal-associated protein. Parkin is identified as a transcriptional repressor of p53 independent of its ubiquitin ligase function. It also has been found that parkin was associated physically with mitochondrial DNA (mtDNA) in proliferating as well as in differentiated SH-SY5Y neuroblastoma cells.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Parkin/PRKN Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | O60260 |

Technical Details

| | |
|------------------|--|
| Immunogen | E.coli-derived human Parkin recombinant protein (Position: I23-K416). Human Parkin shares 82% and 84% amino acid (aa) sequence identity with mouse and rat Parkin, respectively. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Parkin/PRKN Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.