

## Anti-Profilin 1/PFN1 Antibody Picoband® PE Conjugated

Catalog Number: PB9313-PE

### About PFN1

Profilin-1 is a protein that in humans is encoded by the PFN1 gene. The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. PFN1 is mapped to 17p13.2. At high concentrations, profilin prevents the polymerization of actin, whereas it enhances it at low concentrations. This gene is a 140-amino acid protein and major growth regulator of filamentous (F)-actin through its binding of monomeric (G)-actin. It is thought to regulate actin polymerization in response to extracellular signals. Deletion of this gene is associated with Miller-Dieker syndrome.

### Overview

Product Name	Anti-Profilin 1/PFN1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P07737

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Profilin 1, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Profilin 1/PFN1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.