

## Anti-Oct-1/POU2F1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9317-Fluoro550

### About POU2F1

POU domain, class 2, transcription factor 1, also known as OTF1 or OCT1, is a protein that in humans is encoded by the POU2F1 gene. The POU2F1 transcription factor was among the first identified members of the POU transcription factor family. This gene is mapped to 1q24.1 to q24.2. It has been found that the POU dimer from the POU2F1 gene formed on the palindromic OCT factor recognition element, or PORE (ATTTGAAATGCAAAT), could recruit the transcriptional coactivator OBF1, whereas POU dimers formed on the consensus MORE (more PORE) (ATGCATATGCAT) or on MOREs from immunoglobulin heavy chain promoters (AT[G/A][C/A]ATATGCAA) failed to interact. POU domain proteins contain a bipartite DNA-binding domain divided by a flexible linker that enables them to adopt various monomer configurations on DNA.

### Overview

Product Name	Anti-Oct-1/POU2F1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P14859

### Technical Details

Immunogen	E.coli-derived human Oct-1 recombinant protein (Position: S11-Q240). Human Oct-1 shares 96% and 99% amino acid (aa) sequence identity with mouse and rat Oct-1, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Oct-1/POU2F1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.