

Anti-Prolactin/PRL Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9361-Fluoro488

About PRL

Prolactin (PRL), also known as luteotropic hormone or luteotropin, is a protein that in humans is best known for its role in enabling female mammals to produce milk. This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.

Overview

Product Name	Anti-Prolactin/PRL Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01236

Technical Details

Immunogen	E.coli-derived human Prolactin recombinant protein (Position: L29-227C). Human Prolactin shares 60% and 63% amino acid (aa) sequence identity with mouse and rat Prolactin, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

4 Publications Citing This Product

1. PubMed ID: -, Xianglong Meng,Chenzi Lyu,Junnan Ma,Xiaoyan Zhang,Cong Hu,Xiaojuan Su,Chenxu Ning,Wenbin Xie,Shuosheng Zhang,"Metabolomics and Network Pharmacology-Based Investigation into the Mechanisms Underlying the Therapeutic Effect of a New Chinese Traditional Medicine (Cui Nai Ling) on Bromocriptine-Induced Hypogalactia",Evidence-Based Complementary and Alternative Medicine,vol.2021,Article ID 8857449,17 pages,2021.<https://doi.org/10.1155/2021/8857449>
2. PubMed ID: 27827961, Positive Regulation of Decidualization by I-Type Amino Acid Transporter 1 (lat1) in Pregnant Mice
3. PubMed ID: 24133591, MAPK/ERK signal pathway involved expression of COX-2 and VEGF by IL-1? induced in human endometriosis stromal cells in vitro

Visit bosterbio.com/anti-prolactin-picoband-trade-antibody-pb9361-boster.html to see all 4 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prolactin/PRL Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.