

Anti-Prolactin/PRL Antibody Picoband® PE Conjugated

Catalog Number: PB9361-PE

About PRL

Prolactin (PRL), also known as luteotropic hormone or luteotropin, is a protein that in humans is best known for its role in enabling female mammals to produce milk. This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.

Overview

Product Name	Anti-Prolactin/PRL Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01236

Technical Details

Immunogen	E.coli-derived human Prolactin recombinant protein (Position: L29-227C). Human Prolactin shares 60% and 63% amino acid (aa) sequence identity with mouse and rat Prolactin, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

4 Publications Citing This Product

1. PubMed ID: -, Xianglong Meng,Chenzi Lyu,Junnan Ma,Xiaoyan Zhang,Cong Hu,Xiaojuan Su,Chenxu Ning,Wenbin Xie,Shuosheng Zhang,"Metabolomics and Network Pharmacology-Based Investigation into the Mechanisms Underlying the Therapeutic Effect of a New Chinese Traditional Medicine (Cui Nai Ling) on Bromocriptine-Induced Hypogalactia",Evidence-Based Complementary and Alternative Medicine,vol.2021,Article ID 8857449,17 pages,2021.<https://doi.org/10.1155/2021/8857449>

2. PubMed ID: 27827961, Positive Regulation of Decidualization by I-Type Amino Acid Transporter 1 (lat1) in Pregnant Mice

3. PubMed ID: 24133591, MAPK/ERK signal pathway involved expression of COX-2 and VEGF by IL-1? induced in human endometriosis stromal cells in vitro

Visit bosterbio.com/anti-prolactin-picoband-trade-antibody-pb9361-boster.html to see all 4 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prolactin/PRL Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.