

## Anti-BCRP/ABCG2 Antibody Picoband® FITC Conjugated

Catalog Number: PB9364-FITC

### About ABCG2

ABCG2 (Atp-binding cassette, subfamily g, member 2) also known as ABCP, BCRP or MRX, is a protein that in humans is encoded by the ABCG2 gene. The ABCG2 gene encodes a membrane transporter belonging to the ATP-binding cassette (ABC) superfamily of membrane transporters, which are involved in the trafficking of biologic molecules across cell membranes. The ABCG2 protein is also a high capacity transporter for uric acid excretion in the kidney, liver, and gut. The ABCG2 gene is mapped on 4q22.1. In vitro assays of isolated membrane preparations revealed a high-capacity, vanadate-sensitive ATPase activity associated with ABCG2 expression that was stimulated by compounds known to be transported by this protein. Ozvegy et al. (2001) concluded that ABCG2 is likely functioning as a homodimer or homooligomer in this expression system since it is unlikely that putative Sf9 transport partners would be overexpressed at similarly high levels. Abcg2 transports pheophorbide-a, which occurs in various plant-derived foods and food supplements and is highly efficient in limiting its uptake from ingested food. ABCG2 is a major factor in the concentrative transfer of drugs, carcinogens, and dietary toxins to the milk of mice, cows, and humans.

### Overview

Product Name	Anti-BCRP/ABCG2 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal G9
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9UNQ0

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human ABCG2, different from the related mouse sequence by five amino acids, and from the related rat sequence by eight amino acids.
Predicted Reactive Species	Hamster
Cross Reactivity	No cross reactivity with other proteins.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.

Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 7 Publications Citing This Product

1. PubMed ID: 10.3892/ol.2013.1103, Characterization of side population cells isolated from the gastric cancer cell line SGC-7901
2. PubMed ID: 23426065, Li R, Wu X, Wei H, Tian S. Oncol Lett. 2013 Mar;5(3):877-883. Epub 2013 Jan 2. Characterization Of Side Population Cells Isolated From The Gastric Cancer Cell Line Sgc-7901.
3. PubMed ID: 26563263, Reversion of malignant phenotypes of human glioblastoma cells by  $\beta$ -elemene through  $\beta$ -catenin-mediated regulation of stemness-, differentiation- and epithelial-to-mesenchymal transition-related molecules

Visit [bosterbio.com/anti-bcrp-abcg2-picoband-trade-antibody-pb9364-boster.html](http://bosterbio.com/anti-bcrp-abcg2-picoband-trade-antibody-pb9364-boster.html) to see all 7 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BCRP/ABCG2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.