

Anti-AIF/AIFM1 Antibody Picoband® Biotin Conjugated

Catalog Number: PB9366-Biotin

About AIFM1

Apoptosis-inducing factor 1, mitochondrial, also known as AIF or PDCD8 is a protein that in humans is encoded by the AIFM1 gene. AIFM1 gene is mapped to Xq26.1 based on an alignment of the AIFM1 sequence with the genomic sequence. This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells, and it is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it affects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Mutations in this gene cause combined oxidative phosphorylation deficiency 6, which results in a severe mitochondrial encephalomyopathy. A related pseudogene has been identified on chromosome 10.

Overview

Product Name	Anti-AIF/AIFM1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O95831

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human AIF, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

4 Publications Citing This Product

1. PubMed ID: 10.1016/j.pharep.2017.06.005, Icaritin induces apoptosis in acute promyelocytic leukemia by targeting PIM1
2. PubMed ID: 10.1111/j.1478-3231.2008.01757.x, Protection by bicyclol derivatives against acetaminophen-induced acute liver failure in mice and its active mechanism
3. PubMed ID: 27895670, Wnt5a Increases Properties of Lung Cancer Stem Cells and Resistance to Cisplatin through Activation of Wnt5a/PKC Signaling Pathway

Visit bosterbio.com/anti-aif-picoband-trade-antibody-pb9366-boster.html to see all 4 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AIF/AIFM1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.