

Anti-Unrip/STRAP Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9407-Fluoro488

About STRAP

Unrip is also known as STRAP (Serine-threonine kinase receptor-associated protein). It is an enzyme that in humans is encoded by the STRAP gene. It is mapped to 12p12.3. UNRIP is integrated into a complex with UNR or with the SMN complex in vivo in a mutually exclusive manner. It is concluded that UNRIP is the first component of the uridine-rich snRNP assembly machinery that associates with the SMN complex in a compartment-specific way, and it may play a crucial role in the intracellular distribution of the SMN complex.

Overview

Product Name	Anti-Unrip/STRAP Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y3F4

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Unrip, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Unrip/STRAP Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.