

## Anti-Apg7/ATG7 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PB9479-Fluoro594

### About ATG7

Autophagy-related protein 7 is a protein that in humans is encoded by the ATG7 gene. It is mapped to 3p25.3. This gene was identified based on homology to *Pichia pastoris* GSA7 and *Saccharomyces cerevisiae* APG7. In the yeast, the protein appears to be required for fusion of peroxisomal and vacuolar membranes. The protein also shows homology to the ATP-binding and catalytic sites of the E1 ubiquitin activating enzymes. ATG7 is essential for the Apg12 conjugation system that mediates membrane fusion in autophagy. It is found that when nutrients are limited, ATG7 can regulate p53-dependent cell cycle and cell death pathways.

### Overview

Product Name	Anti-Apg7/ATG7 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O95352

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Apg7, different from the related mouse and rat sequences by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

---

1. PubMed ID: 29113167, Yang J, Sheng S, Yang Q, Li L, Qin S, Yu S, Zhang X. Oncol Lett. 2017 Nov;14(5):5333-5339. doi: 10.3892/ol.2017.6857. Epub 2017 Aug 31. Endocan silencing induces programmed cell death in hepatocarcinoma
2. PubMed ID: 29552291, Erbin exerts a protective effect against inflammatory bowel disease by suppressing autophagic cell death
3. PubMed ID: 27383629, Nrf2 signalling and autophagy are involved in diabetes mellitus-induced defects in the development of mouse placenta

Visit [bosterbio.com/anti-apg7-picoband-trade-antibody-pb9479-boster.html](http://bosterbio.com/anti-apg7-picoband-trade-antibody-pb9479-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apg7/ATG7 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.