

## Anti-alpha 1 Sodium Potassium ATPase/ATP1A1 Antibody Picoband® HRP Conjugated

Catalog Number: PB9482-HRP

### About ATP1A1

The protein encoded by this gene belongs to the family of P-type cation transport ATPases, and to the subfamily of Na<sup>+</sup>/K<sup>+</sup> -ATPases. Na<sup>+</sup>/K<sup>+</sup> -ATPase is an integral membrane protein responsible for establishing and maintaining the electrochemical gradients of Na and K ions across the plasma membrane. These gradients are essential for osmoregulation, for sodium-coupled transport of a variety of organic and inorganic molecules, and for electrical excitability of nerve and muscle. This enzyme is composed of two subunits, a large catalytic subunit (alpha) and a smaller glycoprotein subunit (beta). The catalytic subunit of Na<sup>+</sup>/K<sup>+</sup> -ATPase is encoded by multiple genes. This gene encodes an alpha 1 subunit. Multiple transcript variants encoding different isoforms have been found for this gene.

### Overview

Product Name	Anti-alpha 1 Sodium Potassium ATPase/ATP1A1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P05023

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human ATP1A1, different from the related mouse sequence by three amino acids, and from the related rat sequence by four amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

**Suggested Dilutions**

Western blot, Optimal dilutions should be determined by end users.  
Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.  
ELISA, Optimal dilutions should be determined by end users.

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-alpha 1 Sodium Potassium ATPase/ATP1A1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.