

## Anti-Transcription factor Sp4 Antibody Picoband® HRP Conjugated

Catalog Number: PB9558-HRP

### About SP4

SP4 (Sp4 transcription factor), also known as Transcription factor Sp4, MGC130008, MGC130009, SPR-1, HF1B, is a protein that in humans is encoded by the SP4 gene. It is a multigene family of SP proteins. Northern blot analysis that SP4 is abundantly expressed as a 7.5-kb message in brain, but barely detectable in other human cell lines. SP4 is an activator protein like SP1. However, in contrast to SP1, SP4 is not able to act synergistically through adjacent binding sites.

### Overview

Product Name	Anti-Transcription factor Sp4 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q02446

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SP4, different from the related mouse sequence by one amino acid.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor Sp4 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.