

Anti-Cystatin A/CSTA Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9564-Fluoro488

About CSTA

Cystatin A also called Stefin A is a protein that in humans is encoded by the CSTA gene. Cystatin A is a cysteine proteinase inhibitor that belongs to family 1 of the cystatin superfamily. This gene was assigned to human chromosome 3q21.1. This is an intracellular thiol proteinase inhibitor. Cystatin A plays an important role in desmosome-mediated cell-cell adhesion in the lower levels of the epidermis. Stefins have been proposed as prognostic and diagnostic tools for cancer.

Overview

Product Name	Anti-Cystatin A/CSTA Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01040

Technical Details

Immunogen	E.coli-derived human Cystatin A recombinant protein (Position: E8-F98). Human Cystatin A shares 61.5% and 58.2% amino acid (aa) sequence identity with mouse and rat Cystatin A, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cystatin A/CSTA Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.