

Anti-Cystatin A/CSTA Antibody Picoband® HRP Conjugated

Catalog Number: PB9564-HRP

About CSTA

Cystatin A also called Stefin A is a protein that in humans is encoded by the CSTA gene. Cystatin A is a cysteine proteinase inhibitor that belongs to family 1 of the cystatin superfamily. This gene was assigned to human chromosome 3q21.1. This is an intracellular thiol proteinase inhibitor. Cystatin A plays an important role in desmosome-mediated cell-cell adhesion in the lower levels of the epidermis. Stefins have been proposed as prognostic and diagnostic tools for cancer.

Overview

Product Name	Anti-Cystatin A/CSTA Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P01040

Technical Details

Immunogen	E.coli-derived human Cystatin A recombinant protein (Position: E8-F98). Human Cystatin A shares 61.5% and 58.2% amino acid (aa) sequence identity with mouse and rat Cystatin A, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cystatin A/CSTA Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.