

Anti-MCP1/CCL2 Antibody Picoband® FITC Conjugated

Catalog Number: PB9570-FITC

About CCL2

Monocyte chemoattractant protein-1 (MCP-1), a member of the chemokine (chemotactic cytokine) family, is a potent monocyte agonist that is upregulated by oxidized lipids. MCP-1 is also known as CCL2, SCYA2, MCAF. MCAF is a member of family of factors involved in immune and inflammatory responses. The amino acid sequence deduced from the nucleotide sequence reveals the primary structure of the MCAF precursor to be composed of a putative signal peptide sequence of 23 amino acid residues and a mature MCAF sequence of 76 amino acid residues. MCP-1 plays a unique and crucial role in the initiation of atherosclerosis and may provide a new therapeutic target in this disorder. Human MCP-1 is a 8.7KDa non-glycoprotein, consisting of 99 amino acids in precursor form and 76 amino acids in mature form.

Overview

Product Name	Anti-MCP1/CCL2 Antibody Picoband® FITC Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P13500

Technical Details

Immunogen	E. coli-derived human MCP-1 recombinant protein (Position: Q24-T99). Human MCP-1 shares 60.9% and 59.4% amino acid (aa) sequence identity with mouse and rat MCP-1, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

25 Publications Citing This Product

1. PubMed ID: 25210730, Salusin-?, but Not Salusin-?, Promotes Human Umbilical Vein Endothelial Cell Inflammation via the p38 MAPK/JNK-NF- κ B Pathway
2. PubMed ID: 26221302, Comparison of biomarkers in rat renal ischemia-reperfusion injury
3. PubMed ID: 22377751, Liu Z, Wang T, Yang J, Wang S, Yang W, Liu J, Ye Z. Mol Med Rep. 2012 May;5(5):1241-4. Doi: 10.3892/Mmr.2012.813. Epub 2012 Feb 28. Calcium Oxalate Monohydrate Crystals Stimulate Monocyte Chemoattractant Protein-1 And Transforming Grow...

Visit bosterbio.com/anti-mcp-1-picoband-trade-antibody-pb9570-boster.html to see all 25 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MCP1/CCL2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.