

## Anti-Cdc37 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9574-Fluoro550

### About CDC37

Hsp90 co-chaperone Cdc37, also called P50CDC37, is a protein that in humans is encoded by the CDC37 gene. The protein encoded by this gene is highly similar to Cdc 37, a cell division cycle control protein of *Saccharomyces cerevisiae*. This gene is mapped to 19p13.2. And this protein is a molecular chaperone with specific function in cell signal transduction. It has been shown to form complex with Hsp90 and a variety of protein kinases including CDK4, CDK6, SRC, RAF-1, MOK, as well as eIF2 alpha kinases. It is thought to play a critical role in directing Hsp90 to its target kinases.

### Overview

Product Name	Anti-Cdc37 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q16543

### Technical Details

Immunogen	E.coli-derived human Cdc37 recombinant protein (Position: M216-V378). Human Cdc37 shares 95.7% and 95.1% amino acid (aa) sequence identity with mouse and rat Cdc37, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 4 Publications Citing This Product

---

1. PubMed ID: 29404146, Serum Heat Shock Protein Levels and the Relationship of Heat Shock Proteins with Various Parameters in Chronic Obstructive Pulmonary Disease Patients

2. PubMed ID: 23919426, Zhao Z, Miao Y, Pan P, Cheng B, Bai G, Wu H. BMC Complement Altern Med. 2013 Aug 6;13:206. Doi: 10.1186/1472-6882-13-206. Qingfei Xiaoyan Wan Alleviates Asthma Through Multi-Target Network Regulation.

3. PubMed ID: 21346199, Jiang Q, Wang Y, Li T, Shi K, Li Z, Ma Y, Li F, Luo H, Yang Y, Xu C. Mol Biol Cell. 2011 Apr 15;22(8):1167-80. Doi: 10.1091/mbc.E10-10-0860. Epub 2011 Feb 23. Heat Shock Protein 90-Mediated Inactivation Of Nuclear Factor- $\kappa$ B Switches Autophagy To ...

Visit [bosterbio.com/anti-cdc37-picoband-trade-antibody-pb9574-boster.html](http://bosterbio.com/anti-cdc37-picoband-trade-antibody-pb9574-boster.html) to see all 4 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cdc37 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.