

Anti-Apolipoprotein J/CLU Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PB9575-Fluoro594

About CLU

Clusterin (apolipoprotein J) is a 75 - 80 kDa disulfide-linked heterodimeric protein associated with the clearance of cellular debris and apoptosis. In humans, clusterin is encoded by the CLU gene. This protein has several synonyms: dimeric acidic glycoprotein (DAG protein), testosterone repressed prostate message-2 (TRPM-2), sulfated glycoprotein-2 (SGP-2) and complement lysis inhibitor (CLI). Clusterin was mapped to mouse chromosome 14. Clusterin is a ubiquitously expressed molecule thought to influence a variety of processes including cell death. In the brain, it accumulates in dying neurons following seizures and hypoxic-ischemic (H-I) injury. clusterin may be a new therapeutic target to modulate non-caspase-dependent neuronal death following acute brain injury.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Apolipoprotein J/CLU Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P10909 |

Technical Details

| | |
|------------------|--|
| Immunogen | E. coli-derived human Apolipoprotein J recombinant protein (Position: S228-E449). Human Apolipoprotein J shares 76.1 % and 75.2 % amino acid (aa) sequence identity with mouse and rat Apolipoprotein J, respectively. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

2 Publications Citing This Product

1. PubMed ID: 10.1002/2211-5463.12997, Cystatin C regulates the cytotoxicity of infection-induced endothelial-derived beta-amyloid
2. PubMed ID: 33030263, Balczon R, Morrow KA, Leavesley S, Francis CM, Stevens TC, Agwaramgbo E, Williams C, Stevens RP, Langham G, Voth S, Cioffi EA, Weintraub SE, Stevens T. Cystatin C regulates the cytotoxicity of infection-induced endothelial-derived beta-amyloid. FEBS Open Bio

Visit bosterbio.com/anti-apolipoprotein-j-picoband-trade-antibody-pb9575-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein J/CLU Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.