

## Anti-IL8/CXCL8 Antibody Picoband® PE Conjugated

Catalog Number: PB9621-PE

### About CXCL8

Interleukin-8 (IL-8) is a chemokine produced by macrophages and other cell types such as epithelial cells. It is also synthesized by endothelial cells, which store IL-8 in their storage vesicles, the Weibel-Palade bodies. In humans, the interleukin-8 protein is encoded by the IL8 gene. Interleukin-8 (IL8) is a member of the CXC chemokine family. IL-8 is believed to play a role in the pathogenesis of bronchiolitis, a common respiratory tract disease caused by viral infection. This gene and other ten members of the CXC chemokine gene family form a chemokine gene cluster in a region mapped to chromosome 4q. The genes for IL8 have been co-localized on a 335-kb genomic fragment.

### Overview

Product Name	Anti-IL8/CXCL8 Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P10145

### Technical Details

Immunogen	E. coli-derived human IL-8 recombinant protein (Position: A23-S99).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 38 Publications Citing This Product

---

1. PubMed ID: 33550074, Hou J,Lei Z,Cui L,Hou Y,Yang L,An R,Wang Q,Li S,Zhang H,Zhang L.Polystyrene microplastics lead to pyroptosis and apoptosis of ovarian granulosa cells via NLRP3/Caspase-1 signaling pathway in rats.Ecotoxicol Environ Saf.2021 Feb 4;212:112012.doi:10.1016/j.ecotox.2021.112012
2. PubMed ID: 25317156, Neuroprotective effect of pretreatment with ganoderma lucidum in cerebral ischemia/reperfusion injury in rat hippocampus
3. PubMed ID: 25807257, Zhang H, Xia X, Han F, Jiang Q, Rong Y, Song D, Wang Y. Mol Pharm. 2015 May 4;12(5):1648-61. Doi: 10.1021/Acs.Molpharmaceut.5B00069. Epub 2015 Apr 2. Cathelicidin-Bf, A Novel Antimicrobial Peptide From Bungarus Fasciatus, Attenuates Disease In A D...

Visit [bosterbio.com/anti-il-8-picoband-trade-antibody-pb9621-boster.html](http://bosterbio.com/anti-il-8-picoband-trade-antibody-pb9621-boster.html) to see all 38 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL8/CXCL8 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.