

## Anti-Connexin 45/GJA7/GJC1 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9624-Fluoro647

### About GJC1

Gap junction gamma-1 protein (GJC1), also known as gap junction alpha-7 protein (GJA7) or connexin 45 (Cx45), is a protein that in humans is encoded by the GJC1 gene. The International Radiation Hybrid Mapping Consortium mapped the GJA7 gene to chromosome 17q21.31. This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell.

### Overview

Product Name	Anti-Connexin 45/GJA7/GJC1 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P36383

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Connexin 45/GJA7, different from the related mouse and rat sequences by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 27575818, Effect of Novel Gasotransmitter hydrogen sulfide on renal fibrosis and connexins expression in diabetic rats

Visit [bosterbio.com/anti-connexin-45-gja7-picoband-trade-antibody-pb9624-boster.html](https://bosterbio.com/anti-connexin-45-gja7-picoband-trade-antibody-pb9624-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Connexin 45/GJA7/GJC1 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.