

## Anti-Hsp90 beta/HSP90AB1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9636-Fluoro488

### About HSP90AB1

Heat shock protein HSP 90-beta, also called HSP90beta, is a protein that in humans is encoded by the HSP90AB1 gene. It is mapped to chromosome 6p21.1. This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. And this gene is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.

### Overview

Product Name	Anti-Hsp90 beta/HSP90AB1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P08238

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Hsp90 beta, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

---

1. PubMed ID: , Denatured corona proteins mediate the intracellular bioactivities of nanoparticles via the unfolded protein response
2. PubMed ID: -, Biaoxue Rong, Youwen Zhang, Junye Wang et al. Increased STIP1 and Hsp90 Correlate with Progression and Prognosis of Lung Adenocarcinoma, 23 February 2021, PREPRINT (Version 1) available at Research Square [https://doi.org/10.21203/rs.3.rs-226260/v1]

Visit [bosterbio.com/anti-hsp90-beta-picoband-trade-antibody-pb9636-boster.html](https://bosterbio.com/anti-hsp90-beta-picoband-trade-antibody-pb9636-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hsp90 beta/HSP90AB1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.