

## Anti-Hsp70/HSPA1A/HSPA1B Antibody Picoband® PE Conjugated

Catalog Number: PB9638-PE

### About HSPA1A

HSPA1 (heat shock 70kDa protein 1A) also known as HSP70-1, HSPA1A, HSP70-1A, HSP72 or HSP70I, is a protein that in humans is encoded by the HSPA1A gene. This intronless gene encodes a 70kDa heat shock protein which is a member of the heat shock protein 70 family. The HSPA1A gene encodes a predicted 641-amino acid protein. The HSPA1 gene is mapped on 6p21.33. Shimizu et al. (1999) found that peripheral blood mononuclear cells of 18 major depression patients expressed a short HSPA1A transcript that utilized exon 1 rather than exon 2, which is found in the more common HSPA1A transcript. No protein was associated with expression of this short HSPA1A mRNA, possibly due to lack of a TATA box or loss of internal ribosome binding sites. Treatment with BGP-15, a pharmacologic inducer of Hsp72 that can protect against obesity-induced insulin resistance, improved muscular architecture, strength, and contractile function in severely affected diaphragm muscles in mdx dystrophic mice.

### Overview

Product Name	Anti-Hsp70/HSPA1A/HSPA1B Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P0DMV8

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Hsp70, different from the related mouse sequence by five amino acids, and from the related rat sequence by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.

Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 15 Publications Citing This Product

1. PubMed ID: , Denatured corona proteins mediate the intracellular bioactivities of nanoparticles via the unfolded protein response
2. PubMed ID: 31733402, Pan W,Liang J,Tang H,Fang X,Wang F,Ding Y,Huang H,Zhang H.Differentially expressed microRNA profiles in exosomes from vascular smooth muscle cells associated with coronary artery calcification.Int J Biochem Cell Biol.2020 Jan;118:105645.doi:10.1016/j.bioc
3. PubMed ID: 23554704, Effects of minocycline on the expression of NGF and HSP70 and its neuroprotection role following intracerebral hemorrhage in rats

Visit [bosterbio.com/anti-hsp70-picoband-trade-antibody-pb9638-boster.html](http://bosterbio.com/anti-hsp70-picoband-trade-antibody-pb9638-boster.html) to see all 15 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hsp70/HSPA1A/HSPA1B Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.