

Anti-Pea3/ETV4 Antibody Picoband®

Catalog Number: PB9699

About ETV4

ETS translocation variant 4 (ETV4), also known as polyoma enhancer activator 3 (PEA3), or E1AF, is a member of the PEA3 subfamily of Ets transcription factors. It is mapped to 17q21.31. E1AF can activate the promoters of various matrix metalloproteinases, genes whose expression is associated with tumor cell invasion and metastasis, by 10 to 20 fold. Pea3 is detected in cells of epithelial and fibroblastic origin. It is also a transcriptional activator that binds to the enhancer of the adenovirus E1A gene.

Overview


Product Name	Anti-Pea3/ETV4 Antibody Picoband®
Reactive Species	Human, Rat
Description	Boster Bio Anti-Pea3/ETV4 Antibody Picoband® catalog # PB9699. Tested in WB applications. This antibody reacts with Human, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , and 0.05 mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P43268

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Pea3, different from the related mouse sequence by seven amino acids.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human, Rat

Anti-Pea3/ETV4 Antibody Picoband® (PB9699) Images



Western blot analysis of Pea3 using anti-Pea3 antibody (PB9699). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. Lane 1: Rat Lung Tissue Lysate at 50ug, Lane 2: HELA Whole Cell Lysate at 40ug, Lane 3: A549 Whole Cell Lysate at 40ug, Lane 4: SW620 Whole Cell Lysate at 40ug. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Pea3 antigen affinity purified polyclonal antibody (Catalog # PB9699) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for Pea3 at approximately 68 kDa. The expected band size for Pea3 is at 54 kDa.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pea3/ETV4 Antibody

For Research Use Only. Not for use in diagnostic procedures.