

## Anti-MRP8/S100A8 Antibody Picoband® Biotin Conjugated

Catalog Number: PB9742-Biotin

### About S100a8

S100 calcium-binding protein A8 (S100A8) is a protein that in humans is encoded by the S100A8 gene. It is also known as calgranulin A. The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in the inhibition of casein kinase and as a cytokine. Altered expression of this protein is associated with the disease cystic fibrosis.

### Overview

Product Name	Anti-MRP8/S100A8 Antibody Picoband® Biotin Conjugated
Reactive Species	Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P27005

### Technical Details

Immunogen	E.coli-derived mouse S100A8 recombinant protein (Position: P2-E89). Mouse S100A8 shares 58.6% and 80.5% amino acid (aa) sequence identity with human and rat S100A8, respectively.
Cross Reactivity	No cross-reactivity with recombinant human S100A8 is observed.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## 2 Publications Citing This Product

---

1. PubMed ID: , Pretreatment tumor-related leukocytosis misleads positron emission tomography-computed tomography during lymph node staging in gynecological malignancies

2. PubMed ID: 32170086, Mabuchi S,Komura N,Sasano T,Shimura K,Yokoi E,Kozasa K,Kuroda H,Takahashi R,Kawano M,Matsumoto Y,Kato H,Hatazawa J,Kimura T.Pretreatment tumor-related leukocytosis misleads positron emission tomography-computed tomography during lymph node staging in gyne

Visit [bosterbio.com/anti-s100a8-picoband-trade-antibody-pb9742-boster.html](http://bosterbio.com/anti-s100a8-picoband-trade-antibody-pb9742-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MRP8/S100A8 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.