

## Anti-PIGR Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PB9776-Fluoro594

### About PIGR

Polymeric immunoglobulin receptor is a protein that in humans is encoded by the PIGR gene. It is a Fc receptor which facilitates the secretion of the soluble polymeric isoforms of immunoglobulin A and immunoglobulin M. This gene is mapped to 1q31-q41. The encoded poly-Ig receptor binds polymeric immunoglobulin molecules at the basolateral surface of epithelial cells; the complex is then transported across the cell to be secreted at the apical surface. A significant association was found between immunoglobulin A nephropathy and several SNPs in this gene.

### Overview

Product Name	Anti-PIGR Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P01833

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PIGR, different from the related rat sequence by eighteen amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PIGR Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.