

## Anti-Rab5/RAB5A Antibody Picoband® HRP Conjugated

Catalog Number: PB9787-HRP

### About RAB5A

RAB5A (Ras-associated protein RAB5A), also called RAB5, is a protein that in humans is encoded by the RAB5A gene. RAB5 is a rate-limiting component of the machinery regulating the kinetics of membrane traffic in the early endocytic pathway. The RAB5A gene is mapped on 3p24.3. RAB5 is indispensable for a form of receptor tyrosine kinase-induced actin remodeling called circular ruffling. It signals to the actin cytoskeleton through RNTRE, a RAB5-specific GTPase-activating protein (GAP). RAB5 activity on phagosome membranes began to increase on disassembly of the actin coat encapsulating phagosomes. In addition, RAB5 activation is either continuous or repetitive for up to 10 minutes, but it ends before the collapse of engulfed apoptotic cells. Expression of a dominant-negative mutant of RAB5 delayed this collapse of apoptotic thymocytes, showing a role for RAB5 in phagosome maturation.

### Overview

Product Name	Anti-Rab5/RAB5A Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P20339

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Rab5, different from the related mouse and rat sequences by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

**Suggested Dilutions**

Western blot, Optimal dilutions should be determined by end users.  
Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.  
ELISA, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 28491025, Han F, Jiang J, Ding J, Liu H, Xiao B, Shi Y. Front Behav Neurosci. 2017 Apr 26;11:62. doi: 10.3389/fnbeh.2017.00062. eCollection 2017. Change of Rin1 and Stathmin in the Animal Model of Traumatic Stresses

Visit [bosterbio.com/anti-rab5-picoband-trade-antibody-pb9787-boster.html](http://bosterbio.com/anti-rab5-picoband-trade-antibody-pb9787-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rab5/RAB5A Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.