

## Anti-RAB13 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9790-Fluoro647

### About RAB13

Ras-related protein Rab-13 is a protein that in humans is encoded by the RAB13 gene. This gene is a member of the Rab family of small G proteins and plays a role in regulating membrane trafficking between trans-Golgi network (TGN) and recycling endosomes (RE). The encoded protein is involved in the assembly of tight junctions, which are components of the apical junctional complex (AJC) of epithelial cells. The AJC plays a role in forming a barrier between luminal contents and the underlying tissue. Additional functions associated with the protein include endocytic recycling of occludin, regulation of epithelial cell scattering, neuronal regeneration and regulation of neurite outgrowth. Alternately spliced transcript variants have been observed for this gene. A pseudogene associated with this gene is located on chromosome 12.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-RAB13 Antibody Picoband® Fluoro647 Conjugated   |
| Reactive Species     | Human  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | P51153   |

### Technical Details

|                  |   |
|------------------|---|
| Immunogen        | A synthetic peptide corresponding to a sequence at the C-terminus of human RAB13, different from the related mouse and rat sequences by four amino acids. |
| Cross Reactivity | No cross-reactivity with other proteins.  |
| Isotype          | Rabbit IgG  |
| Form             | Liquid  |
| Concentration    | 0.5 mg/mL   |
| Purification     | Immunogen affinity purified.  |
| Conjugate        | Fluoro647<br>Excitation Wavelength: 650 nm<br>Emission Wavelength: 665 nm   |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB13 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.