

Anti-PSD95/DLG4 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9802-Fluoro488

About DLG4

DLG4 discs large homolog 4, also known as PSD95 or SAP-90, is a protein that in humans is encoded by the DLG4 gene. It is a member of the membrane-associated guanylate kinase (MAGUK) family. This gene is mapped to 17p13.1. DLG4 can heteromultimerize with another MAGUK protein, DLG2, and is recruited into NMDA receptor and potassium channel clusters. These two MAGUK proteins may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Overexpression of DLG4 in hippocampal neurons could drive maturation of glutamatergic synapses. DLG4 can orchestrate synaptic development and it has a role in synapse stabilization and plasticity. Ubiquitination of DLG4 through an MDM2-mediated pathway can regulate AMPA receptor surface expression during synaptic plasticity.

Overview

Product Name	Anti-PSD95/DLG4 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P78352

Technical Details

Immunogen	E.coli-derived human PSD95 recombinant protein (Position: H581-L724). Human PSD95 shares 100% amino acid (aa) sequence identity with both mouse and rat PSD95.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm

	Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

2 Publications Citing This Product

1. PubMed ID: 10.1016/j.bbr.2020.113085, Intensive treadmill training promotes cognitive recovery after cerebral ischemia-reperfusion in juvenile rats

2. PubMed ID: -, Ou Qiao,Xinyu Zhang,Yi Zhang et al.Cerebralcare Granule® enhances memantine hydrochloride efficacy in APP/PS1 mice by ameliorating amyloid pathology and cognitive functions,08 April 2021,PREPRINT (Version 1) available at Research Square [https://doi.org/10.21203/rs.3.rs-366097/v1]

Visit bosterbio.com/anti-psd95-picoband-trade-antibody-pb9802-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PSD95/DLG4 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.