

Anti-TAP2 Antibody Picoband® Cy3 Conjugated

Catalog Number: PB9824-Cy3

About TAP2

Transporter, ATP-binding cassette, major histocompatibility complex 2 (TAP2) is a gene in humans that encodes the protein Antigen peptide transporter 2. The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. The gene is assigned to human chromosome 6p21.3. It is located 7 kb telomeric to gene family member ABCB2. The protein encoded by this gene is involved in antigen presentation. And this protein forms a heterodimer with ABCB2 in order to transport peptides from the cytoplasm to the endoplasmic reticulum.

Overview

| | |
|----------------------|--|
| Product Name | Anti-TAP2 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q03519 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human TAP2, different from the related mouse sequence by five amino acids, and from the related rat sequence by six amino acids. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TAP2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.