

Anti-TCP1 alpha Antibody Picoband® PE Conjugated

Catalog Number: PB9826-PE

About TCP1

T-complex protein 1 subunit alpha is a protein that in humans is encoded by the TCP1 gene. The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. In addition, three pseudogenes that appear to be derived from this gene have been found.

Overview

| | |
|----------------------|--|
| Product Name | Anti-TCP1 alpha Antibody Picoband® PE Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P17987 |

Technical Details

| | |
|------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human TCP1, different from the related mouse sequence by one amino acid, and from the related rat sequence by two amino acids. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm |

Suggested Dilutions

Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 16403237, Inhaled nitric oxide alleviates hyperoxia suppressed phosphatidylcholine synthesis in endotoxin-induced injury in mature rat lungs

Visit bosterbio.com/anti-tcp1-picoband-trade-antibody-pb9826-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TCP1 alpha Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.