

Anti-KAP1/TRIM28 Antibody Picoband® APC Conjugated

Catalog Number: PB9835-APC

About TRIM28

Tripartite motif-containing 28 (TRIM28), also known as transcriptional intermediary factor 1beta (TIF1beta) and KAP1 (KRAB-associated protein-1), is a protein that in humans is encoded by the TRIM28 gene. The protein encoded by this gene mediates transcriptional control by interaction with the Kruppel-associated box repression domain found in many transcription factors. The protein localizes to the nucleus and is thought to associate with specific chromatin regions. KAP1 is a ubiquitously expressed protein involved in many critical functions including: transcriptional regulation, cellular differentiation and proliferation, DNA damage repair, viral suppression, and apoptosis. Its functionality is dependent upon post-translational modifications. Phosphorylation of KAP1 acts as a deactivator of the protein in many of its mechanisms while sumoylation acts as an activator.

Overview

Product Name	Anti-KAP1/TRIM28 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q13263

Technical Details

Immunogen	E.coli-derived human KAP1 recombinant protein (Position: A699-P835). Human KAP1 shares 94.9% amino acid (aa) sequence identity with both mouse and rat KAP1.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 33385414, Hao Y, Bai S, Peng J, Hong R, Ding J, Li Z, Guan Y. TRIM27-mediated ubiquitination of PPARgamma promotes glutamate-induced cell apoptosis and inflammation. Exp Cell Res. 2020 Dec 29;112437. doi:10.1016/j.yexcr.2020.112437. Epub ahead of print. PMID:33385414.

Visit bosterbio.com/anti-kap1-picoband-trade-antibody-pb9835-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KAP1/TRIM28 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.