

## Anti-KAP1/TRIM28 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9835-Fluoro488

### About TRIM28

Tripartite motif-containing 28 (TRIM28), also known as transcriptional intermediary factor 1beta (TIF1beta) and KAP1 (KRAB-associated protein-1), is a protein that in humans is encoded by the TRIM28 gene. The protein encoded by this gene mediates transcriptional control by interaction with the Kruppel-associated box repression domain found in many transcription factors. The protein localizes to the nucleus and is thought to associate with specific chromatin regions. KAP1 is a ubiquitously expressed protein involved in many critical functions including: transcriptional regulation, cellular differentiation and proliferation, DNA damage repair, viral suppression, and apoptosis. Its functionality is dependent upon post-translational modifications. Phosphorylation of KAP1 acts as a deactivator of the protein in many of its mechanisms while sumoylation acts as an activator.

### Overview

Product Name	Anti-KAP1/TRIM28 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q13263

### Technical Details

Immunogen	E.coli-derived human KAP1 recombinant protein (Position: A699-P835). Human KAP1 shares 94.9% amino acid (aa) sequence identity with both mouse and rat KAP1.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 33385414, Hao Y,Bai S,Peng J,Hong R,Ding J,Li Z,Guan Y.TRIM27-mediated ubiquitination of PPARgamma promotes glutamate-induced cell apoptosis and inflammation.Exp Cell Res.2020 Dec 29;112437.doi:10.1016/j.yexcr.2020. 112437.Epub ahead of print.PMID:33385414.

Visit [bosterbio.com/anti-kap1-picoband-trade-antibody-pb9835-boster.html](https://bosterbio.com/anti-kap1-picoband-trade-antibody-pb9835-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KAP1/TRIM28 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.