

Anti-UBE2Q2 Antibody Picoband® PE Conjugated

Catalog Number: PB9839-PE

About UBE2Q2

UBE2Q2 is identified as a putative ubiquitin-conjugating enzyme (E2) in a microarray screen for mitotic regulatory proteins. Its gene is mapped to 15q24.2. UBE2Q2 can covalently bind ubiquitin on the active site cysteine within the UBC domain. Inhibition of UBE2Q2 in HeLa cells causes an early mitotic arrest and increases cytotoxicity when the cells are treated with microtubule-inhibiting agents (MIAs). Changes in cell cycle progression and viability are not observed in the absence of MIA treatment, indicating that UBE2Q2 is involved in the response to MIAs rather than performing a more general function in mitosis. Moreover, inhibition of the UBE2Q2 protein causes cells to undergo a prolonged prophase arrest, suggesting that UBE2Q2 normally functions to antagonize an early mitotic checkpoint. Finally, inhibition of UBE2Q2 also sensitizes cells to the cytotoxic effects of MIAs through caspase-mediated apoptosis that is correlated with PARP1 cleavage.

Overview

Product Name	Anti-UBE2Q2 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8WVN8

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human UBE2Q2, different from the related mouse sequence by four amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm

	Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBE2Q2 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.