

## Anti-PGP9.5/UCHL1 Antibody Picoband® FITC Conjugated

Catalog Number: PB9840-FITC

### About UCHL1

UCH-L1, also known as PGP9.5, is a member of a gene family whose products hydrolyze small C-terminal adducts of ubiquitin to generate the ubiquitin monomer. Expression of UCH-L1 is highly specific to neurons and to cells of the diffuse neuroendocrine system and their tumors. It is abundantly present in all neurons (accounts for 1-2% of total brain protein), expressed specifically in neurons and testis/ovary. The catalytic triad of UCH-L1 contains a cysteine at position 90, an aspartate at position 176, and a histidine at position 161 that are responsible for its hydrolase activity.

### Overview

Product Name	Anti-PGP9.5/UCHL1 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P09936

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PGP9.5, different from the related mouse and rat sequences by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

---

1. PubMed ID: 33989617, Zhu WQ,Cai NN,Jiang Y,Yang R,Shi JZ,Zhu CL,Zhang BY,Tang B,Zhang XM.Survivable potential of germ cells after trehalose cryopreservation of bovine testicular tissues.Cryobiology.2021 May

11:S0011-2240(21)00081-X.doi:10.1016/j.cryobiol.2021.05.001.Epub ahead of print.PMID:33989617.

2. PubMed ID: 31305258, Age-associated changes of the intrinsic nervous system in relation with interstitial cells in the pre-weaning goat rumen Yu Liang,1 Imran Tarique,1 Waseem Ail Vistro,1 Yifei Liu,1 Ziyu Wang,2 Abdul haseeb,1 Noor Samad Gandahi,1 Adeela Iqbal,1 Siyi Wang,

3. PubMed ID: 15948235, Nude mice multi-drug resistance model of orthotopic transplantation of liver neoplasm and Tc-99m MIBI SPECT on p-glycoprotein

Visit [bosterbio.com/anti-pgp9-5-picoband-trade-antibody-pb9840-boster.html](https://bosterbio.com/anti-pgp9-5-picoband-trade-antibody-pb9840-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PGP9.5/UCHL1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.