

## Anti-CD14 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9847-Fluoro550

### About CD14

CD14 is a single-copy gene encoding 2 protein forms: a 50- to 55-kD glycosylphosphatidylinositol-anchored membrane protein (mCD14) and a monocyte or liver-derived soluble serum protein (sCD14) that lacks the anchor. This gene is located at bands 5q23-q31. The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. CD14 acts as a co-receptor (along with the Toll-like receptor TLR 4 and MD-2) for the detection of bacterial lipopolysaccharide (LPS). CD14 can bind LPS only in the presence of lipopolysaccharide-binding protein (LBP). Although LPS is considered its main ligand, CD14 also recognizes other pathogen-associated molecular patterns.

### Overview

Product Name	Anti-CD14 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P08571

### Technical Details

Immunogen	E. coli-derived human CD14 recombinant protein (Position: N65-D327). Human CD14 shares 69.8% and 66.8% amino acid (aa) sequence identity with mouse and rat CD14, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 6 Publications Citing This Product

1. PubMed ID: 31930698, Zhou Z,Liu F,Wang L,Zhu B,Chen Y,Yu Y,Wang X.Inflammation has synergistic effect with nicotine in periodontitis by up-regulating the expression of alpha7 nAChR via phosphorylated GSK-3beta.J Cell Mol Med.2020 Feb;24(4):2663-2676.doi:10.1111/jcmm.14986.Epub 2020 Jan 13.PMID:31930698;PMCID:PMC7028870.

2. PubMed ID: 33500729, Shi M,Shen K,Yang B,Zhang P,Lv K,Qi H,Wang Y,Li M,Yuan Q,Zhang Y. An electroporation strategy to synthesize the membrane-coated nanoparticles for enhanced anti-inflammation therapy in bone infection. Theranostics. 2021 Jan 1;11(5):2349-2363.doi:10.7150/th

3. PubMed ID: 23346371, Huang J, Zhu C, Zhang P, Zhu Q, Liu Y, Zhu Z, Wang M, Li W, Yang G, Dong N, Liu J, Chen L, Zhang Y, Yang R, Deng L, Fan J, Wang X, Liu J, Ma B, Fu Q, Wu K. Sci Rep. 2013;3:1114. Doi: 10.1038/Srep01114. Epub 2013 Jan 23. S100+ Cells: A New Neuro-Im...

Visit [bosterbio.com/anti-cd14-picoband-trade-antibody-pb9847-boster.html](https://bosterbio.com/anti-cd14-picoband-trade-antibody-pb9847-boster.html) to see all 6 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD14 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.