

## Anti-ADA Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9914-Fluoro647

### About ADA

Adenosine Deaminase (also known as adenosine aminohydrolase, or ADA) is an enzyme involved in purine metabolism. Primarily, ADA in humans is involved in the development and maintenance of the immune system. However, ADA association has also been observed with epithelial cell differentiation, neurotransmission, and gestation maintenance. It has also been proposed that ADA, in addition to adenosine breakdown, stimulates release of excitatory amino acids and is necessary to the coupling of A1 adenosine receptors and heterotrimeric G proteins. Adenosine deaminase deficiency leads to pulmonary fibrosis, suggesting that chronic exposure to high levels of adenosine can exacerbate inflammation responses rather than suppressing them. It has also been recognized that adenosine deaminase protein and activity is upregulated in mouse hearts that overexpress HIF-1 alpha, which in part explains the attenuated levels of adenosine in HIF-1 alpha expressing hearts during ischemic stress.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-ADA Antibody Picoband® Fluoro647 Conjugated   |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | P00813   |

### Technical Details

|                  |   |
|------------------|---|
| Immunogen        | E. coli-derived human ADA recombinant protein (Position: Q135-L363). Human ADA shares 82.5% and 82.9% amino acid (aa) sequence identity with mouse and rat ADA, respectively. |
| Cross Reactivity | No cross-reactivity with other proteins   |
| Isotype          | Rabbit IgG  |
| Form             | Liquid  |
| Concentration    | 0.5 mg/mL   |
| Purification     | Immunogen affinity purified.  |
| Conjugate        | Fluoro647<br>Excitation Wavelength: 650 nm<br>Emission Wavelength: 665 nm   |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ADA Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.