

Anti-AFG3L2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9915-Fluoro550

About AFG3L2

AFG3L2 is the catalytic subunit of the m-AAA protease, an ATP-dependent proteolytic complex of the mitochondrial inner membrane that degrades misfolded proteins and regulates ribosome assembly. In humans, it is encoded by the AFG3L2 gene. This gene encodes a protein localized in mitochondria and closely related to paraplegin. The paraplegin gene is responsible for an autosomal recessive form of hereditary spastic paraplegia. And this gene is a candidate gene for other hereditary spastic paraplegias or neurodegenerative disorders as well as spastic ataxia-neuropathy syndrome.

Overview

Product Name	Anti-AFG3L2 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y4W6

Technical Details

Immunogen	E. coli-derived human AFG3L2 recombinant protein (Position: R168-D250). Human AFG3L2 shares 100% amino acid (aa) sequence identity with mouse AFG3L2.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AFG3L2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.