

## Anti-Collagen I/COL1A1 Antibody Picoband® FITC Conjugated

Catalog Number: PB9939-FITC

### About Col1a1

Collagen, type I, alpha 1, also known as COL1A1, is a human gene that encodes the major component of type I collagen, the fibrillar collagen found in most connective tissues, including cartilage. This gene is mapped to 17q21.33. And this gene encodes the [pro-alpha1 chains](#) of type I collagen whose triple helix comprises two alpha1 chains and one alpha2 chain. Type I is a fibril-forming collagen found in most connective tissues and is abundant in bone, cornea, dermis and tendon. Mutations in this gene are associated with osteogenesis imperfecta types I-IV, Ehlers-Danlos syndrome type VIIA, Ehlers-Danlos syndrome Classical type, Caffey Disease and idiopathic osteoporosis.

### Overview

Product Name	Anti-Collagen I/COL1A1 Antibody Picoband® FITC Conjugated
Reactive Species	Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P11087

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of mouse Collagen I, different from the related human sequence by three amino acids, and identical to the related rat sequence.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 59 Publications Citing This Product

---

1. PubMed ID: 10.1038/s41598-020-67514-4, Investigating the antifibrotic effect of the antiparasitic drug Praziquantel in in vitro and in vivo preclinical models
2. PubMed ID: 30888602, Wei X, Bao Y, Zhan X, Zhang L, Hao G, Zhou J, Chen Q. MiR-200a ameliorates peritoneal fibrosis and functional deterioration in a rat model of peritoneal dialysis. *Int Urol Nephrol*. 2019 May; 51(5):889-896. doi:10.1007/s11255-019-02122-4. Epub 2019 Mar 19. PMID:30888602; PMCID:PMC6499761.
3. PubMed ID: 31169440, Li X, Bu X, Yan F, Wang F, Wei D, Yuan J, Zheng W, Su J, Yuan J. Deletion of discoidin domain receptor 2 attenuates renal interstitial fibrosis in a murine unilateral ureteral obstruction model. *Ren Fail*. 2019 Nov; 41(1):481-488. doi:10.1080/0886022X.2019.1621759. PMID:31169440; PMCID:PMC6567249.

Visit [bosterbio.com/anti-collagen-i-picoband-trade-antibody-pb9939-boster.html](https://bosterbio.com/anti-collagen-i-picoband-trade-antibody-pb9939-boster.html) to see all 59 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Collagen I/COL1A1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.